

Referencing of NMR spectra, Ξ scale (06/2011)

Several possibilities exist to reference your NMR spectra: Addition of an internal standard (for ^1H and ^{13}C usually TMS), use of an external standard (e.g. for ^{31}P or ^{19}F) with the aid of a capillary tube or, last but not least, by employing the so-called Ξ (Xi) scale [$\Xi_X \equiv 100 \times (\nu_x / \nu_{TMS}^{obs})$], a procedure, which was suggested some years ago by IUPAC and which represents a very useful means of referencing, esp. if heteronuclear spectra are concerned.

In the following, we explain the principle of the Ξ scale approach using the example of an fictitious sample containing ^1H and ^{19}F signals. Note, that the same three step approach can be transferred to any other NMR-active nucleus.

- (1) First of all, a ^1H NMR spectrum needs to be recorded. You need to calibrate the signal of TMS to 0 ppm and then determine the absolute frequency [Hz] of the TMS signal (in our example, this will be done for a 400 MHz ^1H spectrum).
- (2) The frequency value, called $\Theta(\text{TMS})$ or ν_{TMS}^{obs} , read out for the TMS signal (e.g. 400130006.82 Hz for TMS in the ^1H spectrum of our virtual sample) now needs to be multiplied with a reference frequency value, the so-called $\Xi(^{19}\text{F})$, for the heteronucleus in question. In case of ^{19}F , this value is 94.094011 (see table on next page). Finally, the result is divided by 100 MHz, to yield the frequency of the chemical shift of 0 ppm in ^{19}F .

[expressed as an equation: $\nu_F = \nu_{TMS}^{obs} \frac{\Xi_F}{100}$].

- (3) This value (in our example 376.4983726 MHz) can be entered in the NMR software as the value for the reference frequency (in TopSpin, it is the parameter SF , accessible e.g. in the edp menu). Like this, you change the ^{19}F NMR spectrum to a corrected referencing, related to the TMS reference in the ^1H NMR spectrum.

The approach described above is the one favoured by IUPAC, it is highly recommended for use with heteronuclear spectra and easy to apply. On the next page, a table with Ξ values for selected nuclei is provided, which was taken from the original IUPAC publication first published in 2001 in Pure and Applied Chemistry (*Pure Appl. Chem.* **2001**, 73, 1795)

Tabelle 1: Die Spineigenschaften von Spin- $\frac{1}{2}$ -Kernen.^[6]

Isotop ^[6]	natürliche Häufigkeit x [%] ^[6]	magnetisches Moment μ [μ_N] ^[6]	gyromagnetisches Verhältnis γ [10^7 rad s $^{-1}$ T $^{-1}$] ^[6]	Frequenzverhältnis Ξ [%] ^[6]	Referenzverbindung	Probel ^[7]	Lit. für Ξ	D^a	relative Rezeptivität ^[6] D^c
^1H	99.9885	4.837353570	26.7522128	100.000000 ^[6]	Me ₂ Si	CDCl ₃ , $\varphi = 1\%$	–	1.000	5.87×10^3
^2H ^[8]	–	5.159714367	28.5349779	106.663974	[^3H] ₂ Me ₂ Si	[l]	[10]	–	–
^3He	1.37×10^{-4}	–3.685154336	–20.3801587	76.179437	He	Gas	[11]	6.06×10^{-7}	3.56×10^{-3}
^{13}C	1.07	1.216613	6.728284	25.145020	Me ₂ Si	CDCl ₃ , $\varphi = 1\%$	[12, 13]	1.70×10^{-4}	1.00
^{15}N	0.368	–0.49049746	–2.71261804	10.136767	MeNO ₂	in Substanz/CDCl ₃ ^[K]	[9]	3.84×10^{-6}	2.25×10^{-2}
^{19}F	100	4.553333	25.18148	94.094011	CCl ₂ F	[l]	[14]	0.834	4.90×10^3
^{29}Si	4.67	–0.96179	–5.3190	19.867187	Me ₂ Si	CDCl ₃ , $\varphi = 1\%$	[15]	3.67×10^{-4}	2.16
^{31}P	100	1.95999	10.8394	40.480742	H ₃ PO ₄	[l]	[16]	6.65×10^{-2}	3.91×10^2
^{57}Fe	2.119	0.1569636	0.8680624	3.237778	[Fe(CO) ₅]	C ₆ D ₆ ^[l]	[9]	7.24×10^{-7}	4.25×10^{-3}
^{75}Se	7.63	0.92677577	5.1253857	19.071513	Me ₂ Se	in Substanz/C ₆ D ₆ ^[M]	[9]	5.37×10^{-4}	3.15
^{89}Y	100	–0.23801049	–1.3162791	4.900198	Y(NO ₃) ₃	H ₂ O/D ₂ O ^[M]	[9]	1.19×10^{-4}	0.700
^{103}Rh	100	–0.1531	–0.8468	3.186447 ^[6,8]	[Rh(acac) ₃] ^[l]	CDCl ₃ , gesätt.	[18]	3.17×10^{-5}	0.186
(^{107}Ag)	51.839	–0.19689893	–1.0889181	4.047819	AgNO ₃	D ₂ O, gesätt.	[9]	3.50×10^{-5}	0.205
^{109}Ag	48.161	–0.22636279	–1.2518634	4.653533	AgNO ₃	D ₂ O, gesätt.	[9]	4.94×10^{-5}	0.290
(^{111}Cd)	12.80	1.0303729	–5.6983131	21.215480	Me ₂ Cd	in Substanz ^[l]	[19]	1.24×10^{-3}	7.27
^{113}Cd ^[4]	12.22	–1.0778568	–5.9609155	22.193175	Me ₂ Cd	in Substanz ^[l]	[19]	1.35×10^{-3}	7.94
(^{115}Sn)	0.34	–1.5915	–8.8013	32.718749	Me ₂ Sn	in Substanz/C ₆ D ₆ ^[M]	[9]	1.21×10^{-4}	0.711
(^{117}Sn)	7.68	–1.73385	–9.58879	35.632259	Me ₂ Sn	in Substanz/C ₆ D ₆ ^[M]	[9]	3.54×10^{-3}	20.8
^{119}Sn	8.59	–1.81394	–10.0317	37.290632	Me ₂ Sn	in Substanz/C ₆ D ₆ ^[M]	[9]	4.53×10^{-3}	26.6
(^{123}Te)	0.89	–1.276431	–7.059098	26.169742	Me ₂ Te	in Substanz/C ₆ D ₆ ^[M]	[9]	1.64×10^{-4}	0.961
^{125}Te	7.07	–1.5389360	–8.5108404	31.549769	Me ₂ Te	in Substanz/C ₆ D ₆ ^[M]	[9]	2.28×10^{-3}	13.4
^{129}Xe	26.44	–1.347494	–7.452103	27.810186	XeOF ₄	in Substanz ^[l]	[20,21]	5.72×10^{-3}	33.6
^{183}W	14.31	0.20400919	1.1282403	4.166387	Na ₂ WO ₄	D ₂ O, 1 M	[11]	1.07×10^{-5}	6.31×10^{-2}
^{187}Os	1.96	0.1119804	0.6192895	2.282331	OsO ₄	CCl ₄ , 0.98 M	[22]	2.43×10^{-7}	1.43×10^{-3}
^{195}Pt	33.832	1.0557	5.8385	21.496784 ^[6]	[PtCl ₂] ²⁻	D ₂ O, 1.2 M	[9]	3.51×10^{-3}	20.7
^{199}Hg	16.87	0.87621937	4.8457916	17.910822	Me ₂ Hg ^[l]	in Substanz	[11]	1.00×10^{-3}	5.89
(^{209}Tl)	29.524	2.80983305	15.5393338	57.123200 ^[6]	Tl(NO ₃) ₃	[l]	[24]	5.79×10^{-2}	3.40×10^2
^{205}Tl	70.476	2.83747094	15.6921808	57.683838	Tl(NO ₃) ₃	[l]	[25]	0.142	8.36×10^2
^{207}Pb	22.1	1.00906	5.58046	20.920599	Me ₂ Pb	in Substanz/C ₆ D ₆ ^[M]	[9]	2.01×10^{-3}	11.8